



## PRACTICAL WORK N° 5

### THE ACIDIC VINEGAR DEGREE DETERMINATION

#### INTRODUCTION

Vinegar is an aqueous solution with a low content of acetic (ethanoic) acid, which is used as a condiment or food preservative.

#### 2.THE ACIDIC VINEGAR DEGREE (d°) OR (%):

The acidic vinegar degree is indicated on each commercial bottle, and noted (d°) or (%) is: **the mass (g) of pure acetic (ethanoic) acid, contained in 100g of solution (100g vinegar)**. Therefore, the percentage of acidity corresponds to a percentage by mass.

Example: vinegar at 6° or 6% contains 6g of pure acetic (ethanoic) acid per 100g of vinegar solution.

#### Calculation example:

For a volume of 100g of vinegar (Ac: acetic acid)

$$m_{Ac} = n_{Ac} * M_{Ac} = C_{Ac} * V_{Ac} * M_{Ac} \dots (1) \quad (\text{from the formula } n = C * V)$$

We have also :  $d_{Ac} = \rho_{Ac}$  (from the formula  $d_{\text{solution}} = \rho_{\text{solution}} / \rho_{\text{water}}$  ( $\rho_{\text{water}} = 1$ ))

$$= m_{Ac} / V_{Ac}$$

$$V_{Ac} = m_{Ac} / \rho_{Ac} \dots \dots \dots (2)$$

We replace (2) on (1):  $m_{Ac} = m_{Ac} / \rho_{Ac} * C_{Ac} * M_{Ac}$  . So ,  $m_{Ac} = 100 / \rho_{Ac} * C_{Ac} * M_{Ac}$

**The acidic vinegar degree (d°) or (%) =  $100 / \rho_{Ac} * C_{Ac} * M_{Ac}$  (g)**

The determination of the acidic degree of vinegar can be followed by:

#### 2.1 An acid-base titration (as seen on PW3)


A titration reaction is used to determine the concentration of a species in solution ( $C_{Ac} = \dots?$ ) (Ac: Acetic acid)

A first solution of unknown concentration “**titrate solution  $C_{Ac}$** ” is brought into contact with a second solution of known concentration “**titrant solution  $C_B$** ” (B=NaOH). The titrant solution is gradually poured into the titrate solution **drop by drop** until the equivalence point is reached.

**At the equivalence point**  
 $n_{\text{eq.g}}(\text{Ac}) = n_{\text{eq.g}}(\text{NaOH})$

#### 2.2 A pH metric titration:

During a pH-metric titration, the pH of the “**titrate solution  $C_{Ac}$** ” is measured for each volume of “**titrant solution  $C_B$** ” poured. In order to be able to represent regularly the experimental points, the titrant solution must be added milliliter by milliliter.



The pH-metric titration curves ( $pH = f(V_{B\text{added}})$ ) gives the variations in pH as a function of the volume of titrant solution poured show sudden jumps in pH at the equivalence.

The identification of the equivalent volumes, consists of drawing two tangents to the curve  $pH = f(V_{B\text{add}})$ , parallel and placed on either side of the inflection point; then drawing a straight line parallel to these two tangents, equidistant from them. This last straight line intersects the titration curve at the equivalence point, with abscissa  $V_{\text{eq}}$  and ordinate  $pH_{\text{eq}}$ .



### 3.OBJECTIVES:

The main objective of this practical work is :  
Determination of the acidic degrees of commercial vinegar solution.

### 4.PRINCIPAL OF MANIPULATION:

This titration is based on a colorimetric acid-base dosage. And also, on pH metric analysis.

### 5.MATERIAL AND PRODUCTS:

- Burette (25-50 mL)
- Graduated or volumetric pipettes (10 mL)
- Pro-pipette
- Volumetric flasks (100 mL)
- Erlenmeyer flask (100 mL)
- Sodium hydroxide solution (0.1M NaOH)
- Commercial vinegar solution
- pH meter or pH-paper
- Colored indicator
- Beakers

### 6. PREPARATION OF THE DILUTE COMMERCIAL VINEGAR SOLUTION:

1. Take a clean beaker, and pour a small quantity of vinegar solution of unknown concentration  $S_{Ac}$  .
2. Take a volume of 10 ml of the solution  $S_{Ac}$  (use the rinsed pipette).
3. Introduce this test portion into the 100 ml volumetric flask, complete with distilled water, close it with the stopper and shake.

You obtain a solution  $C_{Ac}$  (vinegar solution of unknown concentration).

### 7.TITRATION OF THE DILUTE COMMERCIAL VINEGAR SOLUTION:

#### 7.1Rapid titration (mL by mL) :

- 1- Rinse and Fill the burette with the NaOH solution ( $N_B = 0.1N$ ).
- 2- Take ( $V_{Ac} = 10$  ml) of the dilute vinegar solution and place it in a 100 mL erlenmeyer flask.
- 3- Add 2 drops of the colored indicator.
- 5- Make a rapid titration (mL by mL) to estimate the equivalent volume ( $V_B$ ).
- 6- At the same time, using a pH paper or pH meter, read the pH of the solution  $C_{Ac}$  every mL NaOH until the equivalent volume.
- 7- The solution change color when you added the equivalent volume of titrant solution ( $V_{eq}$ ) . Indicate approximately this volume (by a frame):

$$\dots\dots V_1 \text{ mL} < V_{eq} < \dots\dots V_2 \text{ mL}$$

Table1 : Volumes equivalents

$V_{NaOH}$ titrant (mL)	1	2	3	4	5	6	7	8	9	10
Color solution										

#### 7.2 Precise titration (drop by drop):

- 1- Fill the burette with the NaOH solution ( $N_B = 0.1N$ ).
- 2- Take ( $V_{Ac} = 10$  ml) of the dilute vinegar solution and place it in a 100 mL erlenmeyer flask.
- 3- Add 2 drops of the colored indicator.
- 4- Make precise titration (drop by drop) to estimate the exactly equivalent volume ( $V_{eqB}$ ).
- 5- Indicate the exactly equivalent volume

Table2 : Equivalent volumes :

	1 <sup>st</sup> test	2 <sup>nd</sup> test	$V_{eqB}$ (average)
$V_B$ (mL)			
Color solution			



**8.QUESTIONS:**

- 1-Write down the reaction equation between acetic acid and NaOH.
- 2- Note the exactly equivalent volume  $V_{eqB}$  of (NaOH) added and determine the concentration  $C_{Ac}$  of (dilute vinegar solution).
- 3- Calculate the initial molar concentration  $S_{Ac}$  of the acetic acid (commercial vinegar).
- 4- Calculate the degree of acidity of the vinegar ( $d^\circ$ ) and compare this result with that on the sticker on the commercialized bottle.
- 5- Draw the graph representing the value of  $pH = f(V_B)$  added on the rapid titration (on millimetric paper).
- 6- Determine the equivalent volume  $V_{eqB}$  and the pH at equivalence using the graph.

